

Ambu[®] aScope[™] 4 RhinoLaryngo

Quick Reference Guide



Endoscopic Accessories (Intervention Only)

aScope 4 RhinoLaryngo Intervention: *Accessories up to size 2.0mm may be used through working channel*

Do not use active endoscopic accessories such as laser probes or electro-surgical equipment in conjunction with the aScope 4 Rhinolaryngo system.

Product	Sizes OD/ID	Control Lever Color	Insertion Cord Length
Slim	3.0mm/ -	Purple	300mm
Intervention	5.0mm/ 2.2mm	Green	350mm

Ambu

Image Extend


To extend the live image double click on the live image.
To remove "Image Extend" double click on the image again.

Lubrication

Lubricate the insertion cord starting from the handle and stopping just before the tip with a medical grade lubricant to ensure the lowest possible friction when the aScope 4 RhinoLaryngo is inserted into the patient.


If the camera image of the aScope 4 RhinoLaryngo becomes unclear, clean the tip. When inserting the RhinoLaryngo orally, it is recommended to use a mouthpiece to protect the scope from being damaged.

Instillation of Fluids (Intervention Only)



Insert a syringe into the working channel port at the top of the aScope 4 Rhinolaryngo Intervention (If using a luer lock adapter insert the included introducer on top of the working port. Make sure to insert the port completely. Failure to do so may result in fluid spilling from the working channel port).


Use the syringe to instill fluid. Make sure suction is not applied during this process, as this will



Flush the channel with 2 ml of air to ensure all fluid has left the channel.


CAUTION: Remove the introducer and/or syringe before applying suction. Failure to do so may result in a reduction of suction and/or spillage of fluids.

Insertion of Endoscopic Accessories (Intervention Only)



Insert the endoscopic accessory into the working channel port and advance it carefully through the working channel until it can be seen on the live image on the monitor. The enclosed introducer can be used to facilitate the insertion of soft accessories.

Suction (Intervention Only)



Secure the tubing properly to the suction connector before suction is applied.

CONSIDERATION: In some cases, the suction tubing may need to be cut to obtain a tight seal around the suction connector. Failure to secure the suction connector properly may result in a reduction of suction.



Ambu Inc.
6721 Columbia Gateway Drive
Suite 200
Columbia, MD 21046
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com

Questions? Contact your Ambu representative
at 800-262-8462 or visit www.ambuUSA.com

US Rx Only