Instructions for Use

Cuff Pressure Gauge

CE 0123

EN - English

Instructions for Use 4 - 6

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
Einsteinstrasse 1 | 72172 Sulz a.N. | Germany

Intended Use

The Cuff Pressure Gauge is used for inflation, deflation and pressure control of Tracheal Tubes with high volume low pressure cuffs and Ambu® Aura Laryngeal Mask.

ATTENTION

- Read and follow the instructions for use carefully before using the product.
- The product must only be used by trained medical personnel.
- Before use, the products must be inspected visually for damages (cracks, breakage, etc.) and a functional check must be performed (see section "Functional check").
- The cuff pressure gauge is not suitable for MRI!
- In case that the cuff pressure gauge is used for monitoring and remains on the inflation line, the pressure must be checked regularly every 15 minutes by the user. Pressure increase can be adjusted by pressing the red release valve and pressure decrease can be adjusted by squeezing the inflation bulb. This cuff pressure gauge is a manual system and can not hold the cuff pressure constantly over a longer period.
- When disconnecting the cuff pressure gauge from the inflation line, it has to be disconnected together with the connecting tubing as otherwise the system is open.
- A mechanical stop is integrated in the housing to prevent that the release valve gets stuck. The release valve can get stuck by excessive strong activation. In this case, the device must not be used further and must be returned to Ambu for repair.
- The connecting tube is for single use and must not be reprocessed. Reprocessing influences the function of the connecting tube. Any reuse involves the potential risk for cross contamination.

Functional Check

Test the Cuff Pressure Gauge as follows:

- 1a) Close Luer nozzle with the finger.
- 1b) Inflate with the inflation bulb to 40 cmH₂O; the value must be constant for 2 - 3 sec; in case of pressure decrease, the device must be returned to Ambu for repair.

If a connecting tubing is used, the system consisting of Cuff Pressure Gauge and connecting tubing has to be checked as follows (additionally to 1):

- 2a) Attach the connecting tubing to the Luer nozzle of the Cuff Pressure Gauge.
- 2b) Seal the end of the connecting tubing with the finger.
- 2c) Inflate with the hand inflation bulb to 40 cmH₂O; the value must be constant for 2 - 3 sec; in case of a pressure decrease, the connecting tubing is leaking and must be replaced. Check again with the new connecting tubing (repeat steps 2a - 2c).

During connection and disconnection of the connecting tube, do not bend the luer fitting of the device, to prevent damage. Furthermore, make sure that the connecting tubing is held on the connector (not at the tubing).

Instruction for use with Tracheal Tubes with high-volume low-pressure cuff

Before use check the low pressure cuff for leaks. Prior to intubation or extubation withdraw all air from the cuff via a syringe. Connect the Cuff Pressure Gauge (eventually by using the connecting tubing) to the inflation line of the Tracheal Tube. Inflate the cuff to a pressure of 60 - 90 cmH₂O. This ensures that the cuff is in close contact with the tracheal wall. Release the pressure immediately by pressing the red release valve until the pointer reaches the green area "Tracheal tube". We recommend the intra cuff pressure to be above 22 cmH₂O due to the risk of aspiration or pneumonia and below 32 cmH₂O due to the risk of ischemia of the tracheal mucous membrane (see publication of J. Stauffer in Respiratory Care, July 1999).

Check and monitor the cuff pressure continuously. Adjust pressure increase via the red release valve and pressure decrease via the inflation bulb.

Use with Ambu® Aura Laryngeal Mask

Functional testing as described in the directions for use must be carried out before using the Ambu Aura Laryngeal Mask. Insert the Ambu Aura Laryngeal Mask according to directions as described in the directions for use. Attach the Ambu Cuff Pressure Gauge and its connecting tubing to the inflation line of the cuff. Without holding the tube, inflate the cuff by pumping with just enough air to obtain a seal, equivalent to intracuff pressures of approximately 60 cmH₂O. The pressure in the cuff can be checked constantly by connecting the Ambu Cuff Pressure Gauge. In case of a pressure increase press the red release valve to adjust to the desired pressure. The pressure should not be higher than 60 cmH₂O.

Check and monitor the cuff pressure continuously. Adjust pressure increase via the red release valve and pressure decrease via the inflation bulb.

Extubation

Deflate the cuffs completely.

Wipe disinfection

Wipe disinfection using conventional surface disinfectants based on alcohol (e.g. Incides® N from Ecolab Deutschland GmbH) is suitable for disinfection. When choosing a product for disinfection, a disinfectant with appropriate ranges of action must be used: bactericidal, yeasticidal, tuberculocidal and virucidal. After the wipe disinfection, the devices must be inspected for visible contamination. If necessary, repeat the wipe disinfection. After the wipe disinfection, the devices must be checked in accordance with the section "Functional check".

CAUTION

- Do not immerse the product in liquids.
- The product cannot be automatic or manual cleaned respectively sterilised.
- Do not clean the glass (plastic material) with aggressive agents.

Metrological Inspection

The accuracy ± 2 cmH₂O of the Cuff Pressure Gauges must be checked every 24 months. If required meas-

uring equipment is not available, the metrological control can be carried out by Ambu. In this case the Cuff Pressure Gauge has to be sent to Ambu.

If there is an indication that the Cuff Pressure Gauge cannot keep the accuracy of measurement within the 24 months, or the metrological characteristics of the Cuff Pressure Gauge have been affected, a metrological control has to be carried out immediately. This has to be carried out only by a person who meets the requirements for metrological controls. The results of the control, as well as the measured values and measurement procedures has to be documented.

After successful control, the Cuff Pressure Gauge has to be marked with a sign which indicates (clearly and traceable) the year of the next inspection and the person who carried out the metrological control.

In order to check the values, the Luer nozzle of the Cuff Pressure Gauge has to be connected to a calibrated pressure gauge (measurement inaccuracy max. ± 0.6 cmH₂O). Afterwards inflate via the inflation bulb to 30, 60 and 90 cmH₂O and record the values. The values must be within the indicated tolerance. If this is not the case, the device must be returned to Ambu for repair. If the values are within the tolerance, the device can still be used.

Note for conversion of cmH₂O to SI-Unit Pascal:

$$1 \text{ cmH}_2\text{O} = 0,98067 \text{ hPa}$$

Service

Prior to return of medical products for complaint / repair, they must be cleaned and disinfected to exclude any risk for Ambu's staff. For safety reasons Ambu preserves the right to refuse contaminated products.

Storage and transport conditions

The shelf life of the Cuff Pressure Gauge is 8 years.

The shelf life of the connecting tube is 5 years.

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

Maintenance

The medical product may only be serviced by persons who have the technical knowledge, qualification and the required means to properly carry out the service.

After service of the medical product, the essential constructive and functional characteristics for safety and functionality have to be checked, insofar as they can be influenced by the service procedure.

To ease repair of the medical product, return it together with a detailed description of the defect.

Disposal

Used or damaged products must be disposed in accordance with the applicable national and international regulations.

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