

Ambu® VivaSight™ 2 DLT

THE SECOND GENERATION OF AMBU DOUBLE LUMEN TUBES FOR ONE LUNG VENTILATION

VivaSight 2 DLT

Increase patient safety with continuous visualization

VivaSight 2 DLT can increase patient safety because it enables you to follow tube placement and detect and correct tube malpositioning.

Improve your procedure flow

Continuous visualization can positively affect patient management and procedure workflow because it:


- Allows faster intubation time than conventional DLTs.
- Facilitates correct tube placement.
- Significantly reduces the use of a bronchoscope.
- Allows immediate detection and correction of tube dislodgement.

A cost-effective alternative for OLV

According to a health-economic study based on a randomized controlled trial, VivaSight DLT offers a cost-effective alternative for OLV procedures and demonstrates potential savings of \$47.65 per procedure.

The study compared the costs and effects of a conventional DLT used in combination with a reusable bronchoscope to video-enabled DLT alone.

VivaSight 2 DLT can be connected with Ambu's full HD monitor, aView 2 Advance, for high quality imaging and Dual View.



KEY BENEFITS

Increases patient safety

Improves procedure workflow

Cost-effective alternative for One Lung Ventilation

Significantly reduces the use of a bronchoscope

Compatible with Ambu® aView™ 2 Advance with Dual View mode

Ambu

SPECIFICATIONS

Ambu® VivaSight™ 2 DLT	35 Fr Left	37 Fr Left	39 Fr Left	41 Fr Left
Size	35 Fr	37 Fr	39 Fr	41 Fr
ID - Effective Inside Diameter Bronchial (mm)	4.4	4.6	4.8	5.0
ID - Effective Inside Diameter Tracheal (mm)	4.4	4.6	4.8	5.0
OD (mm- bronchial)	10.5	11.0	11.5	12.0
Length (±10 mm)	330	330	330	330

VivaSight 2 DLT Kit Includes	Quantity
VivaSight 2 DLT with Y-connector	1
Ambu® VivaSight™ Suction Catheter	2
VivaSight 2 Adapter Cable	1

General	
Usage	Single-use only
Pack	Supplied individually packaged, ETO-sterilized
Shelf life	36 months
Storage	Store at temperatures between 10 oC - 25 oC in a dry, cool and dark place

Materials	
Tube	Medical-grade PVC; additional injection port for rinsing camera optics
Cuff	Medical-grade PVC, high volume, low pressure
Tube hardness	85 (Shore A)
This product is not made with phthalates	

Camera	
Image sensor: CMOS	

Light Source	
LEDs: 2 white	


ORDERING SPECIFICATIONS

Ambu® VivaSight™ 2 DLT Kits	Order number	Quantity/box
Ambu® VivaSight 2 DLT Kit 35 Fr - LEFT	412351000	5
Ambu® VivaSight 2 DLT Kit 37 Fr - LEFT	412371000	5
Ambu® VivaSight 2 DLT Kit 39 Fr - LEFT	412391000	5
Ambu® VivaSight 2 DLT Kit 41 Fr - LEFT	412411000	5

Monitors	Order number	Quantity/box
Ambu® aView™ 2 Advance Monitor	405011000	1



Ambu® VivaSight™ 2 DLT



Ambu® VivaSight™ 2 DLT connected to Ambu® aView™ 2 Advance Monitor



Ambu, Inc.
 6721 Columbia Gateway Drive, Suite 200
 Columbia, Maryland 21046
 Tel. 800 262 8462
 Fax 800 262 8673
 ambuUSA.com

CE 2797 US: Rx only