

# Ambu® AuraStraight™ Silicone Single-Use Laryngeal Mask

## DATASHEET

The AuraStraight Silicone is intended for use as an alternative to the face mask for achieving and maintaining control of the airway during routine and emergency anesthetic procedures in fasted patients.

AuraStraight Silicone is a single-use laryngeal mask. The absence of epiglottis bars ensures easier tracheal access.

The mask cuff is designed to conform to the contours of the hypopharynx and is available in seven sizes, ranging from neonate through adult.



## KEY BENEFITS

Single-use

Clearly marked guidelines for easy insertion

Cuff inflation volume, and weight recommendations listed on the airway tube

Available in seven sizes from neonate to adult

PVC and phthalate-free material

# Ambu

## SPECIFICATIONS

	Infant		Pediatric		Adult		
Mask size	#1	#1½	#2	#2½	#3	#4	#5
Patient weight	<5 kg	5-10 kg	10-20 kg	20-30 kg	30-50 kg	50-70 kg	70-100 kg
Maximum cuff inflation volume	4 ml	7 ml	10 ml	14 ml	20 ml	30 ml	40 ml
Maximum intra-cuff pressure	60 cm H <sub>2</sub> O						
Airway connector	15 mm male (ISO)						
Min. I.D. tube	5.2 mm	6.6 mm	7.0 mm	8.5 mm	9.8mm	9.8 mm	11.5 mm
Max. O.D. tube	8.0 mm	9.8 mm	10.6 mm	12.7 mm	14.8 mm	14.8 mm	16.5 mm
Inflation valve	Luer cone (DIN 13090 F6)						

NOTE: Appropriate storage temperature -40°C to 60°C



AuraStraight Silicone laryngeal mask

## MATERIAL

Product	Material
Airway tube	Silicone
Cuff	Silicone
Airway connector	Polycarbonate
Pilot tube	Silicone
Pilot balloon with inflation valve	Silicone/PP/SS

Packaging	Quantity	Material
Pouch	1	Tyvek®/Polyethylene Terephthalate (PET) Polyethylene (PE)
Box	10	Cardboard

## ORDERING SPECIFICATIONS

Item No.	Size	Type
342 100 000	1	Single-Use Straight Laryngeal Mask
342 150 000	1 ½	Single-Use Straight Laryngeal Mask
342 200 000	2	Single-Use Straight Laryngeal Mask
342 250 000	2 ½	Single-Use Straight Laryngeal Mask
342 300 000	3	Single-Use Straight Laryngeal Mask
342 400 000	4	Single-Use Straight Laryngeal Mask
342 500 000	5	Single-Use Straight Laryngeal Mask

Questions? Contact your Ambu representative at **800-262-8462** or visit **ambuUSA.com**

**Ambu** FOREVER FORWARD

Ambu, Inc.  
6721 Columbia Gateway Drive  
Suite 200  
Columbia, MD 21046  
Tel. 800 262 8462  
Fax 800 262 8673  
ambuUSA.com

US Rx only

LBL-008917 - V02 - 02/2024 - Ambu USA. Technical data may be modified without further notice.