

VENTILATE YOUR PATIENTS WITH CONFIDENCE

The SPUR II from Ambu combines optimal tidal volume delivery* with capnography to help ensure safe & reliable manual resuscitation. Utilizing its unique M-port and a gas sampling line, the SPUR II facilitates sidestream EtCO₂ monitoring by conveniently providing caregivers with all the tools necessary to properly ventilate patients in an emergency situation.

Ambu® SPUR® II with EtCO₂



KEY BENEFITS

Produces clean and consistent waveforms

Unique filtration system helps protect the monitor from moisture

Direct M-port connection reduces dead space

Eliminates need for separate sample line and adapters resulting in SKU reduction

Reflective connector compatible with most capnography monitors

Delivers tidal volumes which meet AHA standards for cardiopulmonary resuscitation*

Product	Material
520211000E	Ambu Adult SPUR II Bag Reservoir, Medium Mask, Gas Sampling Line
520211000BE	Ambu Adult SPUR II Bag Reservoir, Medium Mask, Gas Sampling Line, Individually Boxed
520211001E	Ambu Adult SPUR II Bag Reservoir, Medium Mask, PEEP Valve, Gas Sampling Line
520211010E	Ambu Adult SPUR II Bag Reservoir, Medium Mask, 10' O ₂ Tubing, Gas Sampling Line
520211061E	Ambu Adult SPUR II Bag Reservoir, Medium Mask, Expiratory Filter, PEEP Valve, Gas Sampling Line
520214000E	Ambu Adult SPUR II Bag Reservoir, Medium & Large Masks, Gas Sampling Line
520611000E	Ambu Adult SPUR II Bag Reservoir, Medium Mask, Gas Sampling Line
520611048E	Ambu Adult SPUR II Bag Reservoir, Medium Mask, Expiratory Filter, Gas Sampling Line
521611000E	Ambu Adult SPUR II Tube Reservoir, Medium Mask, Gas Sampling Line
523211000E	Ambu Adult SPUR II Bag Reservoir, Pop-off Valve, Medium Mask, Manometer Port, Gas Sampling Line
523211001E	Ambu Adult SPUR II Bag Reservoir, Pop-off Valve, Medium Mask, Manometer Port, PEEP Valve, Gas Sampling Line
523611030E	Ambu Adult SPUR II Bag Reservoir, Pop-off Valve, Medium Mask, Manometer, Gas Sampling Line
523611031E	Ambu Adult SPUR II Bag Reservoir, Pop-off Valve, Medium Mask, Manometer, PEEP Valve, Gas Sampling Line
523611051E	Ambu Adult SPUR II Bag Reservoir, Pop-off Valve, Medium Mask, Expiratory Filter, Manometer, PEEP Valve, Gas Sampling Line

*AHA guidelines state that tidal volumes of approximately 500-600 mL (6-7 mL per kg) are sufficient for adult resuscitation www.ahajournals.org/doi/full/10.1161/CIRCULATIONAHA.110.970939

Ambu

Ambu Inc.
6230 Old Dobbin Lane
Columbia, MD 21045
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com

US: Rx only