


# Ambu® Disposable Pressure Manometer



## KEY BENEFITS

The slim design makes it easy to monitor pressure with minimum obstruction of the view of the patient's face.

The single patient use concept offers minimized risk of cross infection.

This product does not contain any natural rubber latex or PVC.

The manometer has a color-coded scale to give enhanced visual support during ventilation. The color coding scale is only for visual support. The correct ventilation pressure must be determined by the medical professional.

7-Tesla MR conditional

## DATASHEET

The Ambu Disposable Pressure Manometer for single patient use is intended to be used for monitoring the patient's airway pressure.

The Ambu Disposable Pressure Manometer is suitable for use with Ambu resuscitators or other resuscitators, hyperinflation bags, CPAP masks or circuits.



Ambu Disposable Pressure Manometer connected to Ambu SPUR II Resuscitator

# Ambu

## SPECIFICATIONS

### Size

Length	55 mm (2.17")
Diameter	22 mm (0.87")
Weight	6.4 g
Connector size	ID 4.0 mm (0.16")

### Materials

Connector	SEBS
Spring	CuSn
Housing	SB
Piston	PP
Membrane	Silicone
Operating temperature range	-18°C a 50°C
Storage temperature range	-40°C a 60°C

### Measuring

Measuring Limits	0-60 cm H <sub>2</sub> O (hPa)
Accuracy	+/- 2 cm H <sub>2</sub> O (hPa) a 7.5, 10, 15, 20, 30 H <sub>2</sub> O (hPa) +/- 3 cm H <sub>2</sub> O (hPa) a 40 H <sub>2</sub> O (hPa) +/- 5 cm H <sub>2</sub> O (hPa) a 60 H <sub>2</sub> O (hPa)



Ambu Disposable Pressure Manometer connected to Ambu SPUR II Pediatric

## ORDERING SPECIFICATIONS

### Item No.

322 004 000	Ambu Disposable Pressure Manometer US
-------------	---------------------------------------

Questions? Contact your Ambu representative at **800-262-8462** or visit **[www.ambuUSA.com](http://www.ambuUSA.com)**

**Ambu** FOREVER FORWARD

Ambu Inc.  
6721 Columbia Gateway Drive  
Suite 200  
Columbia, MD 21046  
Tel. 800 262 8462  
Fax 800 262 8673  
[ambuUSA.com](http://ambuUSA.com)

US Rx Only