

Kit i.o. Access



Key features

- Training of i.o. access procedures
- Realistic anatomical landmarks
- Easy to use
- Reusable injection pad
- May be used with all devices for i.o.

Kit i.o. Access

The Ambu Kit i.o. Access allows realistic training of intraosseous infusion practice.

Promotes training scenarios where the trainees can learn to identify the palpable landmarks including the patella, tibia, and tibial tuberosity.

The Ambu Kit i.o. Access is compatible with all commercially available i.o. devices. The Ambu Kit i.o. is applied using realistic pressure.

Designed for realistic training.

Simple maintenance and economical follow up costs

Well thought down to the smallest details, Ambu provides a realistic i.o. training solution which is easily maintainable. The Kit i.o. Access allows smart replacement of the consumables like bone, foam and skin.

Compatible to devices for establishing i.o. access

It can be used with all commercially available devices for establishing intraosseous access.

Easy upgrade for existing training manikins

Ambu provides a smart upgrade for your current Ambu training manikin to include the i.o. access possibility.

Materials

Leg	PVC (Polyvinylchlorid) hard
Bone	HDPE (High density Polyethylen)
Injection pad	IXPE Polyethylen
Skin	PVC (Polyvinylchlorid) soft

Order information

Item no	Description
284 000 520	Kit i.o. Access incl. 1 bone, 5 injection pads, 5 skins




Kit i.o. Access



Kit i.o. Access



Realistic anatomic landmarks



Ambu Inc.

6230 Old Dobbin Lane
Columbia, MD 21045
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com