

BRONCHOSCOPY LIKE YOU'VE NEVER SEEN IT BEFORE

Ambu® aScope™ 5 Broncho combines the excellent maneuverability and next-level imaging needed in the bronchoscopy suite with the sterility and efficiency of the single-use concept. The result? Consistently high performance, sterility, and the flexibility you need for a complex and demanding workflow. It is compatible with the full-HD resolution Ambu® aBox™ 2 displaying and processing unit, which can easily be connected to external monitors.



aScope 5 Broncho delivers consistently high performance and excellent maneuverability, supported by a new rotational function and a smart functional design.

EXCELLENT MANEUVERABILITY

Access difficult-to-reach areas of the bronchial tree:

The combination of the 195°/195° up-and-down bending angles, distal tip angulation range, and insertion tube flexibility helps you access difficult-to-reach areas in the bronchial tree. You will notice the excellent maneuverability, especially when working with an endotherapy instrument.

120° left/right rotation enables you to adjust the position of the distal end to ease navigation and insertion of endotherapy instruments.





Ambu® aScope™ 5 Broncho

NEXT-LEVEL IMAGING

High-resolution camera chip and full-HD displaying and processing unit (12.8" screen): aScope 5 Broncho and aBox 2 work together to deliver exceptional imaging.

Excellent image quality:

Advanced image processing, 120° field-of-view and 3-100 mm depth-of-field HD system with advanced adaptive auto control give you the bright and clear image needed for the procedure.

Easy to document images and videos

At the touch of a button on the bronchoscope, you can zoom in and out, capture images and video, and review them with patients or colleagues post-procedure. Documentation is smooth and efficient thanks to compatibility with PACS servers and electronic medical record software.

Contact your Ambu Representative to schedule a demo at **800-262-8462** or visit **www.ambuUSA.com**

Ambu FOREVER
FORWARD

Ambu, Inc.
6230 Old Dobbin Lane, Suite 250
Columbia, MD 21045
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com