

**STERILE.  
SINGLE-USE.  
PROBLEM SOLVED.**

Ambu® aScope™ Duodeno

**Ambu**



## STERILE, SINGLE-USE

Patient cross-contamination has been a problem. Sterile straight from the pack, Ambu aScope Duodeno is a single-use duodenoscope for ERCP procedures. It gives you a new scope with reliable performance every time, and there's no risk of patient cross-contamination.

When it comes to sterility, cost-effectiveness and reliable performance, the Ambu aScope Duodeno single-use endoscope for ERCP is a game-changer.



### STERILE STRAIGHT FROM THE PACK

aScope Duodeno is sterile, eliminating the risk of scope to patient cross-contamination.

### NO MORE REPROCESSING OR REPAIR

With single-use, there is no need for post-procedure cleaning or repair, thus optimizing workflow and productivity.

### PERFORMANCE YOU CAN TRUST

The familiar design ensures a seamless transition from conventional scopes and delivers the performance you expect.



# Explore the new aScope Duodeno solution

Produced by the world's leading single-use endoscope supplier

Ambu launched the world's first single-use flexible bronchoscope in 2009, used in more than 6,000 hospitals worldwide. aScope Duodeno builds on a decade of know-how placing a reliable, sterile and cost-effective solution in your hands.

## CONTROL AND STABILITY IN YOUR HANDS

The ergonomic, ripple-surface handle and lightweight design of the aScope Duodeno is comfortable to use, giving you the same control and stability as your current reusable duodenoscopes.



## INNOVATION-DRIVEN DESIGN

With high-definition imagery and flexible bending angles, you are able to view the mucosa in detail and navigate the gastrointestinal tract efficiently. Our precise elevator function allows you to access the biliary system with confidence.



## A SIMPLE SOLUTION

aScope Duodeno is part of a complete endoscope solution consisting of the single-use duodenoscope and our cost-effective and multi-functional Ambu® aBox™ Duodeno processor.





Item code	Product	Quantity
<b>aScope Duodeno</b>		
482001000 Canada Only 482101000 US Only	aScope Duodeno	6 scopes
<b>aBox Duodeno</b>		
485001000US	aBox Duodeno	1 processor

## Ambu takes environmental issues seriously

Danish companies are renowned for prioritizing sustainability. Ambu aScope Duodeno eliminates the need to use gases and liquid chemical germicides byproducts that can impact air and water quality. Inspired by this, Ambu has set ambitious sustainable design goals and actively seeks out potential partnerships with service providers to find ways to help you dispose of our products in the most cost-effective, risk-free and environmentally-friendly way possible.

## A history of breakthrough ideas

Ambu has been bringing the solutions of the future to life since 1937. Today, millions of patients and healthcare professionals worldwide depend on the efficiency, safety and performance of our single-use endoscopy, anaesthesia, and patient monitoring & diagnostics solutions.

Headquartered near Copenhagen, Denmark, Ambu employs approximately 4,000 people in Europe, North America and the Asia Pacific region.

Contact your local territory manager to schedule a demo  
[ambuUSA.com/single-use-duodenoscope](https://www.ambuUSA.com/single-use-duodenoscope)



Ambu USA  
6721 Columbia Gateway Drive, Suite 200  
Columbia, Maryland 21046  
Tel. 800 262 8462  
Fax 800 262 8673  
[ambuUSA.com](https://www.ambuUSA.com)

CE 2797 US: Rx only