

DATASHEET

Sterile straight from the pack, Ambu aScope Duodeno is a single-use endoscope with side viewing optics, deflectable tip and an elevator to control the position of compatible endoscopic accessories.

With aScope Duodeno, you have a new scope with consistent performance every time, and there is no risk of patient cross-contamination.

What's more, there's no need for costly reprocessing or repair.



Distal tip showing elevator in use with a sphincterotome

Ambu® aScope™ Duodeno

**STERILE.
SINGLE-USE.
PROBLEM SOLVED.**



KEY BENEFITS

Sterile, single-use



No risk of patient cross-contamination

Cost transparency

Eliminates the risk and cost of reprocessing

Familiar design and function to that of conventional duodenoscopes

Ambu



SPECIFICATIONS - AMBU® ASCOPE™ DUODENO

Optical System	
Direction of view	Retro side viewing 6°
Field of view	130°
Insertion section	
Distal end outer diameter	φ13.7 mm
Insertion tube outer diameter	φ11.3 mm
Working length	1240 mm
Instrument working channel	
Inner diameter	φ4.2 mm (12.6 FR)
Bending section	
Angulation range	Up 120° / Down 90° Left 90° / Right 110°
Illumination method	Dual LEDs
Dimensions	
Weight	700 g

ORDERING SPECIFICATIONS*

Reference Number	Description	Quantity
485001000	aBox Duodeno	1
482001000 Canada Only	aScope Duodeno	Box of 6
482001000/4 Canada Only	aScope Duodeno	Box of 4**
482101000 US Only	aScope Duodeno	Box of 6
482101000/4 US Only	aScope Duodeno	Box of 4**

*Some configurations may not be available worldwide

**The outer box contains four inner boxes, each containing one scope.

Questions? Contact your Ambu Representative
at 800-262-8462 or visit www.ambuUSA.com

Ambu FOREVER FORWARD

Ambu, Inc.
6721 Columbia Gateway Drive, Suite 200
Columbia, Maryland 21046
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com

CE 2797 US: Rx only