


DATASHEET

Ambu Subdermal Hooked Needle and Ambu Subdermal Curved Needle are stainless steel electrodes, which are designed to stay in place during usage. The electrodes have color coded lead wires which aid the identification of electrodes. The products come with a 1.5 mm touch proof DIN 42802 connectors and can be used with most multi-channel nerve monitors.

The hooked and curved needle electrodes are developed to be used for monitoring, recording and stimulation of bio and nerve potential signals.

Ambu® Subdermal Hooked Needle

Ambu® Subdermal Curved Needle



KEY FEATURES

- Designed to stay better in place during usage
- Color coded lead wires aid the identification of electrodes
- Needle electrode made of stainless steel
- 1.5 mm touch proof DIN 42802 connector

RECOMMENDED APPLICATIONS

- Intraoperative Monitoring (IOM)
- Electroencephalogram (EEG)
- Electromyography (EMG)
- Evoked Potentials (EP)

Ambu

SPECIFICATIONS

Electrode

Connector type	1.5 mm touch proof DIN 42802
Sterilization method	Ethylene oxide (EO)

Environment

Electrode, lead wire and packaging are not made with natural rubber latex
Packaging is PVC-free


MATERIALS

Electrode


Needle	Stainless Steel
Lead wire	PVC insulated Tin Plated Copper Wire
Shrink tube	Polyolefin
Protection tube	Silicone

Packaging

Pouches, transparent layer	Polyester/Polyethylene
Pouches, paper layer	Paper
Boxes	Cardboard



Subdermal Hooked Needle



Subdermal Curved Needle

ORDERING SPECIFICATIONS

Item no.	Lead wire length		Lead wire color	Shelf life in months (unopened pouches)	Packaging		
	cm	inch			Units/pouch	Units/inner box	Units/outer box
Subdermal Hooked Needle							
74912-120/24	120	47"	○●●●●●	36	1	24	1152
74912-200/24	200	79"		36	1	24	1152
74912-250/24	250	98"		36	1	24	1152
Subdermal Curved Needle							
75021-150/25	150	59"	○●●●●	36	1	25	1200
75021-200/25	200	79"		36	1	25	1200

Questions? Contact your Ambu representative at **800-262-8462** or visit **ambuUSA.com**



Ambu Inc.
6721 Columbia Gateway Drive
Suite 200
Columbia, MD 21046
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com



US Rx Only